

Keep the "O"s in the Nose!



It's all about the babies®

Sticky Whiskers®

Being secure has never been so close



Sticky Whiskers®:

- *help keep cannula in place directly under infant's nose with soft, low-profile Velcro®*
- *attaches with gentle hydrocolloid adhesive*
- *easy to change placement of cannula*
- *one size fits all*
- *can be used for securing OG/NG tubes*



BMS
Beever's Manufacturing & Supply

Telephone: 503.472.9055
Toll Free: 800.818.4025
Fax: 503.434.6303
www.beevers.net
info@beevers.net



When securement matters



Sticky Whiskers® :

- attaches with gentle, non-irritating hydrocolloid adhesive
- secures cannula in place with soft, low-profile Velcro®
- allows for easy repositioning of cannula

Request a free sample now!

Size & Description

One size fits all babies

Product #

SW102

Beevers Manufacturing & Supply

800.818.4025 • 503.434.6303 (fax) • www.beevers.net • info@beevers.net